Package 'RPA'

March 26, 2013

Type Package

Version 1.14.0 Date 2012-09-08 Author Leo Lahti < leo.lahti@iki.fi> Maintainer Leo Lahti < leo.lahti@iki.fi> Depends R (>= 2.15.1), affy, affydata, methods, parallel biocViews GeneExpression, Microarray, Preprocessing, QualityControl Description Probabilistic analysis of probe reliability and differential gene expression on short oligonucleotide arrays. Lahti et al. "Probabilistic Analysis of Probe Reliability in Differential Gene Expression Studies with Short Oligonucleotide Arrays", TCBB/IEEE, 2011.	
Author Leo Lahti <leo.lahti@iki.fi> Maintainer Leo Lahti <leo.lahti@iki.fi> Depends R (>= 2.15.1), affy, affydata, methods, parallel biocViews GeneExpression, Microarray, Preprocessing, QualityControl Description Probabilistic analysis of probe reliability and differential gene expression on short oligonucleotide arrays. Lahti et al. "Probabilistic Analysis of Probe Reliability in Differential Gene Expression Studies with</leo.lahti@iki.fi></leo.lahti@iki.fi>	
Maintainer Leo Lahti <leo.lahti@iki.fi> Depends R (>= 2.15.1), affy, affydata, methods, parallel biocViews GeneExpression, Microarray, Preprocessing, QualityControl Description Probabilistic analysis of probe reliability and differential gene expression on short oligonucleotide arrays. Lahti et al. "Probabilistic Analysis of Probe Reliability in Differential Gene Expression Studies with</leo.lahti@iki.fi>	
Depends R (>= 2.15.1), affy, affydata, methods, parallel biocViews GeneExpression, Microarray, Preprocessing, QualityControl Description Probabilistic analysis of probe reliability and differential gene expression on short oligonucleotide arrays. Lahti et al. "Probabilistic Analysis of Probe Reliability in Differential Gene Expression Studies with	
biocViews GeneExpression, Microarray, Preprocessing, QualityControl Description Probabilistic analysis of probe reliability and differential gene expression on short oligonucleotide arrays. Lahti et al. "Probabilistic Analysis of Probe Reliability in Differential Gene Expression Studies with	
Description Probabilistic analysis of probe reliability and differential gene expression on short oligonucleotide arrays. Lahti et al. "Probabilistic Analysis of Probe Reliability in Differential Gene Expression Studies with	
differential gene expression on short oligonucleotide arrays. Lahti et al. "Probabilistic Analysis of Probe Reliability in Differential Gene Expression Studies with	
http://doi.ieeecomputersociety.org/10.1109/TCBB.2009.38	
License FreeBSD	
LazyLoad yes	
R topics documented:	
RPA-package estimate.affinities estimate.hyperparameters get.batches get.probe.noise.estimates hyperparameter.update online.quantile plot probe.performance rpa rpa-class rpa.fit rpa.fit-class	2 3 4 5 6 7 9 10 11 12 14 15 16

2 RPA-package

RPA	-package	RP	4: _I	oro	ba	bil	list	ic	an	al	ysi	s	of p	ore	ob	e r	eli	iał	oili	ty	ar	ıd	ge	ne	e.	хр	re.	ssi	on	
Index																														31
	summarize.batches			•		•			•		•		٠		•		•	•	•		•		•	•			•			29
	sample.probeset																													28
	s2.update																													27
	rpa2eset																													26
	RPA.preprocess																													25
	RPA.pointestimate																													23
	rpa.plot																													22
	rpa.online																													20
	rpa.list-class									•			•	•	•		•	•			•		٠	•		•	•	•		19

Description

RPA estimates probe-specific variances and differential gene expression using probe-level observations of differential gene expression.

Details

Package: RPA
Type: Package
Version: 1.7.32
Date: 2011-03-12
License: FreeBSD
LazyLoad: yes

RPA.pointestimate is used to calculate probe reliability and differential expression estimates: 'rpa.results <- RPA.pointestimate(affybatch)'. The other functions are provided for users who wish to investigate the details of the algorithm more closely.

Author(s)

Leo Lahti <leo.lahti@iki.fi>

References

Probabilistic Analysis of Probe Reliability in Differential Gene Expression Studies with Short Oligonucleotide Arrays. Lahti et al., TCBB/IEEE 2011.

Examples

```
## Not run:

## Load example data set (Dilution affybatch).

## This is a toy example with a small example dataset

## for probe reliability analysis (4 arrays).

## For practical applications, a larger sample size is
```

estimate.affinities 3

```
## recommended.

#require(affy)
#require(affydata)
#data(Dilution)

## Compute RPA for whole data set
## ... slow, not executed here
#rpa.results <- RPA.pointestimate(Dilution)
```

estimate.affinities

Probe affinity estimation

Description

Estimates probe-specific affinity parameters.

Usage

```
estimate.affinities(dat, mu)
```

Arguments

dat Input data set: probes x samples.

mu Estimated expression signal from RPA model.

Details

To estimate means in the original data domain let us assume that each probe-level observation x is of the following form: x = d + mu + noise, where x and d are vectors over samples, mu is a scalar (vector with identical elements) noise is Gaussian with zero mean and probe-specific variance parameters tau2 Then the parameter mu will indicate how much probe-level observation deviates from the estimated signal shape d. This deviation is further decomposed as mu = mu.real + mu.probe, where mu.real describes the 'real' signal level, common for all probes mu.probe describes probe affinity effect Let us now assume that $mu.probe \sim N(0, sigma.probe)$. This encodes the assumption that in general the affinity effect of each probe tends to be close to zero. Then we just calculate ML estimates of mu.real and mu.probe based on particular assumptions. Note that this part of the algorithm has not been defined in full probabilistic terms yet, just calculating the point estimates.

Note that while tau2 in RPA measures stochastic noise, and NOT the affinity effect, we use it here as a heuristic solution to weigh the probes according to how much they contribute to the overall signal shape. Intuitively, probes that have little effect on the signal shape (i.e. are very noisy and likely to be contaminated by many unrelated signals) should also contribute less to the absolute signal estimate. If no other prior information is available, using stochastic parameters tau2 to determine probe weights is likely to work better than simple averaging of the probes without weights. Also in this case the probe affinities sum close to zero but there is some flexibility, and more noisy probes can be downweighted.

Value

A vector with probe-specific affinities.

Note

Affinity estimation is not part of the original RPA procedure in TCBB/IEEE 2011 paper. It is added here since estimates of the absolute levels are often needed in microarray applications. Note that affinity parameters are unidentifiable in the model if no prior assumptions are given. We assume that affinity effects are zero on average, but allow flexibility through probe-specific weights.

Author(s)

Leo Lahti <leo.lahti@iki.fi>

References

See citation("RPA").

See Also

rpa.fit

Examples

```
\#\# mu <- estimate.affinities(dat, mu)
```

estimate.hyperparameters

Hyperparameter estimation

Description

Hyperparameter estimation

Usage

estimate. hyperparameters (sets = NULL, priors = list(alpha = 2, beta = 1), batches, cdf = NULL, quantile. basis

Arguments

sets Probesets to handle. All probesets by default.

priors User-defined priors

batches Data batches for online learning cdf CDF probeset definition file quantile.basis Basis for quantile normalization bg.method Background correction method

epsilon Convergence parameter

load.batches Logical. Load preprocessed data whose identifiers are picked from names(batches).

Assuming that the same batch list (batches) was used to create the files in on-

line.quantiles function.

save. hyperparameter. batches

Save hyperparameters for each batch into files using the identifiers with batch

name with -hyper.RData suffix.

get.batches 5

mc.cores Number of cores for parallel computation

verbose Print progress information

 $normalization. \\ method$

Normalization method

save.batches.dir Specify the output directory for temporary batch saves.

unique.run.identifier

Define identifier for this run for naming the temporary batch files. By default, a

random id is generated.

Details

Sweeps through the batches. Updates hyperparameters at each batch.

Value

alpha Hyperparameter alpha (same for all probesets)

betas Hyperparameter beta (probe-specific) variances Probe-specific variances (beta/alpha)

Author(s)

Leo Lahti <leo.lahti@iki.fi>

References

See citation("RPA").

Examples

#

get.batches	Split data into batches

Description

Splits data into batches

Usage

```
get.batches(items, batch.size, shuffle = FALSE)
```

Arguments

items A vector of items to be splitted into batches.

batch.size Batch size. The last batch may contain less elements than the other batches

which have batch.size elements each.

shuffle Split the elements randomly in the batches.

Value

A list. Each element corresponds to one batch and contains a vector listing the elements in that batch.

Author(s)

Leo Lahti <leo.lahti@iki.fi>

References

See citation("RPA").

See Also

RPA.preprocess

Examples

#

get.probe.noise.estimates

Fetch probe-level noise estimates from an rpa object

Description

Provides estimates of probe-level noise, given by the RPA algorithm.

Usage

get.probe.noise.estimates(rpa.res, sets = NULL, normalization = NULL, verbose = FALSE)

Arguments

rpa.res An rpa object. sets Probesets to check.

normalization Optional normalization for probe noise estimates.

The higher the value, the higher the probe-level noise. By default, probe-level variances of the RPA model are returned. Other options include:

"withinset.weights": The relative weight of a probe within probeset is determined by the relative noise of the probe with respect to the other probes in the same probeset. This option returns the inverse of probe-specific weights within each probeset. This can be used to normalize probe-level weights to improve comparability across probesets.

"withinset.relative": The detected probe-level noise can be coupled with overall signal levels of the probeset. This option provides an estimate of probe-wise standard deviation normalized by the standard deviation of the probeset-level signal d.

"withinset.categorical": In some applications it can be sufficient to investigate the relative order of the probes, ignoring the parameter estimates. This option indexes the probes according to their reliability within each probeset. Probes with higher indices are more noisy. verbose Print progress information during computation.

Details

The normalization options are included to improve comparability across probesets. The higher the variance, the more noisy the probe. Inverse of the variance, can be used to quantitate probe reliability. Note that the relative weight of a probe within probeset is determined by the relative noise of the probe with respect to the other probes in the same probeset. Comparison of probespecific variances across probesets may benefit from normalization of this effect. Therefore optional normalizations for probe noise estimation are provided.

Value

A list. Each element corresponds to one probeset (of the input object). The element lists noise estimates for each probe within the probeset.

Author(s)

Leo Lahti <leo.lahti@iki.fi>

References

Probabilistic Analysis of Probe Reliability in Differential Gene Expression Studies with Short Oligonucleotide Arrays. Lahti et al., TCBB/IEEE 2011. See http://bioconductor.org/packages/release/bioc/html/RPA.htm

See Also

RPA.pointestimate

Examples

```
## Load example data set require(affydata) data(Dilution) ## Compute RPA rpa.results <- RPA.pointestimate(Dilution, set = "1000_at") noise <- get.probe.noise.estimates(rpa.results)
```

hyperparameter.update Update hyperparameters.

Description

Update shape (alpha) and scale (beta) parameters of the inverse gamma distribution.

Usage

```
hyperparameter.update(dat, alpha, beta, th = 0.01)
```

Arguments

dat A probes x samples matrix (probeset).

alpha, beta Shape and scale parameters of inverse gamma density for the probe variances.

th Convergence threshold.

Details

Shape update: alpha <- alpha + T/2; Scale update: beta <- alpha * s2 where s2 is the updated variance for each probe (the mode of variances is given by beta/alpha). The variances (s2) are updated by EM type algorithm, see s2.update.

Value

A list with elements alpha, beta (corresponding to the shape and scale parameters of inverse gamma distribution, respectively).

Author(s)

```
Leo Lahti <leo.lahti@iki.fi>
```

References

```
See citation("RPA").
```

See Also

s2.update, rpa.online

Examples

```
# Generate and fit toydata, learn hyperparameters
set.seed(11122)
P < -11 \# number of probes
N <- 5000 \# number of arrays
real <- sample.probeset(P = P, n = N, shape = 3, scale = 1, mu.real = 4)
dat <- real\$dat \# probes x samples
# Set priors
alpha < -1e-2
beta <- rep(1e-2, P)
# Operate in batches
step <- 1000
for (ni in seq(1, N, step)) {
 batch <- ni:(ni+step-1)
 hp <- \ hyperparameter.update(dat[,batch], \ alpha, \ beta, \ th = 1e-2)
 alpha <- hp\alpha
 beta <- hp\$beta
\# Final variance estimate
s2 <- beta/alpha
\# Compare real and estimated variances
```

online.quantile 9

```
plot(sqrt(real\$variance), sqrt(s2), main = cor(sqrt(real\$variance), sqrt(s2))); abline(0,1)
```

online.quantile Quantile normalization tools for online preprocessing.

Description

Estimate quantiles for quantile normalization based on subset of the data (random, or specified by the user).

Usage

```
online.quantile(abatch, n) quorm.basis.online(cel.files, bg.method = "rma", cdf = NULL, save.batches = FALSE, batch.size = 2, verbose = TRUE, save.batches.dir = ".", unique.run.identifier = NULL)
```

Arguments

abatch	AffyBatch						
n	Numeric: number of random samples to use to define quantile basis. Vector: specify samples to be used in quantile basis calculation.						
cel.files	List or vector of CEL files to process. Can also be a list of batches, each batch containing a vector of CEL files						
bg.method	Background correction method						
cdf	Optional. CDF file for alternative probeset definitions						
save.batches	Logical. If TRUE, save results for each batch into file with batchname.RData where 'batchname' is picked from the batch identifiers in the batch list (as given by get.batches function).						
batch.size	Batch size						
verbose	Print progress information						
save.batches.dir	Specify the output directory for temporary batch saves.						
unique.run.identifier							
	Define identifier for this run for naming the temporary batch files. By default, a						

Details

"online.quantile": Ordinary quantile normalization is exhaustively memory-consuming in alrge data sets. Then the quantiles can be calculated based on subset of the data to allow efficient normalization. This function can also be used to investigate effect of subset size to convergence of the quantile estimates; "qnorm.basis.online": sweeps through the data in batches to calculate the basis for quantile normalization (average over sorted profiles).

random id is generated.

Value

"online.quantile": AffyBatch; "qnorm.basis.online": a vector containing the basis for quantile normalization.

10 plot

Author(s)

```
Leo Lahti <leo.lahti@iki.fi>
```

References

```
See citation("RPA").
```

See Also

RPA.preprocess

Examples

#

plot

Plot RPA results and probe-level data for a specific probeset.

Description

Plots the preprocessed probe-level observations, estimated probeset-level signal, and probe-specific variances. It is also possible to highlight individual probes and external summary measures.

Usage

```
\#\# S3 method for class 'rpa' plot(x, y = NULL, set, highlight.probes = NULL, pcol = "darkgrey", mucol = "black", ecol = "red", cex.lab = 1.5, cex.axis = 1, cex.names = 1, cex.main = 1, external.sig
```

Arguments

x,y Instances of the 'rpa class; y is optional and never used. Provided for consis-

tency.

set Probeset to plot.

highlight.probes

Optionally highlight some of the probes (with dashed line)

pcol Color for probe signal visualization.

mucol Color for summary estimate.
ecol Color for external signal.
cex.lab, cex.axis, cex.names, cex.main

Axis adjustment parameters.

external.signal Plot external signal on the probeset. For instance, an alternative summary esti-

mate from another preprocessing methods.

main Title for plot.

plots Type of plot. See help(rpa.plot).
... Other arguments to be passed.

probe.performance 11

Value

Used for side-effects. Returns probes x samples matrix of probe-level data plotted on the image.

Author(s)

```
Leo Lahti <leo.lahti@iki.fi>
```

References

Probabilistic Analysis of Probe Reliability in Differential Gene Expression Studies with Short Oligonucleotide Arrays. Lahti et al., TCBB/IEEE. See http://bioconductor.org/packages/release/bioc/html/RPA.html

See Also

RPA.pointestimate

Examples

```
# Not run:

## Load example data set

#require(affydata)

#data(Dilution)

## Compute RPA for specific probesets only

#set <- "1000_at"

#rpa.results <- RPA.pointestimate(Dilution, set)

## Visualize the results for one of the probe sets

#plot(rpa.results, set)
```

probe.performance

probe.performance

Description

Provide a table of probe-level parameter estimates (affinity and stochastic noise) for rpa class object.

Usage

```
probe.performance(rpa.object, abatch = NULL, sets = NULL)
```

Arguments

rpa.object Object of the rpa class (output from function RPA.pointestimate)
abatch Optional: Affybatch for the rpa object, if not provided in the rpa project.

sets Optional: Specify the inspected probesets.

Value

Data frame of probe-level parameter estimates

12 rpa

Author(s)

Leo Lahti <leo.lahti@iki.fi>

References

See citation("RPA")

Examples

library(affydata) data(Dilution); rpa.results <- RPA.pointestimate(Dilution); tab <- probe.parameters(rpa.results)

rpa

RPA for preprocessing.

Description

Returns an expressionSet object preprocessed with RPA. If 'cind' is not specified, uses the first array of affybatch as the reference.

Usage

rpa(abatch=NULL, sets=NULL, priors=NULL, epsilon=1e-2, cind=1, tau2.method="robust", d.method=1e-2, cind=1, tau2.method=1e-2, cind=1, tau3.method=1e-2, cind=1, tau3.method=1e-2, cind=1, tau3.method=1e-2, cind=1e-2, ci

Arguments

abatch An AffyBatch object.

sets Probesets for which RPA will be computed. Default: all probe sets.

priors An 'rpa.priors' object. Can be used to set user-specified priors for the model

parameters. Not used tau2.method = "var". The prior parameters alpha and beta are prior parameters for inverse Gamma distribution of probe-specific variances. Noninformative prior is obtained with alpha, beta -> 0. Not used with tau2.method 'var'. Scalar alpha and beta specify an identical inverse Gamma prior for all probes, which regularizes the solution. Can be also specified as

lists, each element corresponding to one probeset.

epsilon Convergence tolerance. The iteration is deemed converged when the change in

all parameters is < epsilon.

cind Specify reference array for computing probe-level differential expression. De-

fault: cind = 1. Note that if exclude.reference.array = TRUE the expression value for the reference array (cind) will be excluded in the output. Note that all values of the reference array are 0 since they indicate the differential expression

of the reference array against itself.

tau2.method Optimization method for tau2 (probe-specific variances). This parameter is de-

noted by tau^2 in the vignette and manuscript.

"robust": (default) update tau2 by posterior mean, regularized by informative priors that are identical for all probes (user-specified by setting scalar values for alpha, beta). This regularizes the solution, and avoids overfitting where a single probe obtains infinite reliability. This is a potential problem in the other tau2 update methods with non-informative variance priors. The default values alpha

= 2; beta = 1 are used if alpha and beta are not specified.

rpa 13

"mode": update tau2 with posterior mean "mean": update tau2 with posterior mean

"var": update tau2 with variance around d. Applies the fact that tau2 cost func-

tion converges to variance with large sample sizes.

d.method Method to optimize d.

"fast": (default) weighted mean over the probes, weighted by probe variances

The solution converges to this with large sample size.

"basic": optimization scheme to find a mode used in Lahti et al. TCBB/IEEE;

relatively slow; this is the preferred method with small sample sizes.

verbose Print progress information during computation.

bg.method Specify background correction method. Default: "rma". See bgcorrect.methods()

for other options.

normalization.method

Specify quantile normalization method. Default: "pmonly". See normalize.methods(Dilution)

for other options.

cdf Specify an alternative CDF environment. Default: none.

cel.files List of CEL files to preprocess. cel.path Path to CEL file directory.

Details

RPA preprocessing function. Gives an estimate of the probeset-level mean parameter d of the RPA model, and returns these in an expressionSet object.

Value

An instance of the 'expressionSet' class.

Note

The choices tau2.method = "robust" and d.method = "fast" are recommended. With small sample size and informative prior, d.method = "basic" may be preferable. For very large expression data collections, see rpa.online function.

Author(s)

Leo Lahti <leo.lahti@iki.fi>

References

See citation("RPA").

See Also

rpa.online, RPA.pointestimate, set.priors, AffyBatch, ExpressionSet, estimate.affinities, rpa.fit

Examples

```
# Not run:
## Load example data set
#require(affydata)
#data(Dilution)
## eset <- rpa(Dilution)
```

14 rpa-class

rpa-class

Class "rpa"

Description

Class for the RPA package.

Objects from the Class

Returned by RPA.pointestimate function. Objects can be created by calls of the form new ("rpa", ...).

Slots

Contains the following information for analyzed probesets and data:

Object of class "list"

- .DdtaA matrix of probesets x arrays. Specifies the estimated 'true' underlying gene expression signal over the arrays for each investigated probeset.
- **mu.real** Numeric. Used for technical purposes. Indicates the difference between original signal levels and the signal levels in the differential signal domain used to learn the probe-specific variance parameters. See model description in vignette for details.
- **tau2** A list. Each element corresponds to a probeset, and contains a vector that gives the estimated variance (noise level) for each probe in that probeset.
- **affinity** A list. Each element corresponds to a probeset, and contains a vector that gives the estimated affinity effect for each probe in that probeset.
- **cind** Specifies which of the arrays in abatch was used as a reference for computing probe-level differential expression.
- sets A character vector listing the investigated probesets.
- **cdf** Alternative CDF that was used in the analysis.
- data Preprocessed probe-level data on which the model was fitted.
- abatch The associated affybatch object.

Extends

```
Class "list", from data part. Class "vector", by class "list", distance 2.
```

Methods

```
[ signature(x = "rpa"): ...
[[ signature(x = "rpa"): ...
show signature(x = "rpa"): ...
```

Author(s)

Leo Lahti <leo.lahti@iki.fi>

References

Probabilistic Analysis of Probe Reliability in Differential Gene Expression Studies with Short Oligonucleotide Arrays. Lahti et al., TCBB/IEEE. http://www.roihu.info/publications/preprints/TCBB10.pdf

rpa.fit

Examples

```
showClass("rpa")
```

rpa.fit

rpa.fit

Description

Fits the RPA model, including estimation of probe-specific affinity parameters.

Usage

```
rpa.fit(dat, cind = 1, epsilon = 1e-2, alpha = NULL, beta = NULL, tau2.method = "robust", d.method = "fast")
```

Arguments

dat Original data: probes x samples.

cind Index of reference array.

epsilon Convergence tolerance. The iteration is deemed converged when the change in

all parameters is < epsilon.

alpha, beta Priors for inverse Gamma distribution of probe-specific variances. Noninforma-

tive prior is obtained with alpha, beta -> 0. Not used with tau2.method 'var'. Scalar alpha and beta are specify equal inverse Gamma prior for all probes to

regularize the solution. The defaults depend on the method.

tau2.method Optimization method for tau2 (probe-specific variances).

"robust": (default) update tau2 by posterior mean, regularized by informative priors that are identical for all probes (user-specified by setting scalar values for alpha, beta). This regularizes the solution, and avoids overfitting where a single probe obtains infinite reliability. This is a potential problem in the other tau2 update methods with non-informative variance priors. The default values alpha

= 2; beta = 1 are used if alpha and beta are not specified.

"mode": update tau2 with posterior mean

"mean": update tau2 with posterior mean

"var": update tau2 with variance around d. Applies the fact that tau2 cost func-

tion converges to variance with large sample sizes.

d.method Method to optimize d.

"fast": (default) weighted mean over the probes, weighted by probe variances

The solution converges to this with large sample size.

"basic": optimization scheme to find a mode used in Lahti et al. TCBB/IEEE;

relatively slow; preferred with small sample size.

Details

First learns a point estimate for the RPA model in terms of differential expression values w.r.t. reference sample. After this, probe affinities are estimated by comparing original data and differential expression shape, and setting prior assumptions concerning probe affinities.

16 rpa.fit-class

Value

mu Fitted signal in original data: mu.real + d
mu.real Shifting parameter of the reference sample

tau2 Probe-specific stochastic noise

affinity Probe-specific affinities
data Probeset data matrix
alpha, beta prior parameters

Note

Affinity estimation is not part of the original RPA procedure in TCBB/IEEE 2011 paper. It is added here since estimates of the absolute levels are often needed in microarray applications. Note that affinity parameters are unidentifiable in the model if no prior assumptions are given. We assume that affinity effects are zero on average, but allow some flexibility through probe-specific weights.

Author(s)

Leo Lahti <leo.lahti@iki.fi>

References

Probabilistic Analysis of Probe Reliability in Differential Gene Expression Studies with Short Oligonucleotide Arrays. Lahti et al., TCBB/IEEE 2011. See http://www.roihu.info/publications/preprints/TCBB10.pdf

See Also

rpa, RPA.pointestimate, estimate.affinities

Examples

res <- rpa.fit(dat, cind, epsilon, alpha, beta, tau2.method, d.method, affinity.method)

rpa.fit-class *Class "rpa.fit"*

Description

Class for the RPA package.

Objects from the Class

Returned by rpa.fit function. Objects can be created by calls of the form new("rpa.fit", ...).

RPA.iteration 17

Slots

Contains the following information for analyzed probesets and data:

Object of class "list"

.Datu: A matrix of probes x samples. Gives the summary estimate corresponding to the 'true' underlying signal over the samples.

mu.real Numeric. Used for technical purposes. Indicates the difference between original signal levels and the signal levels in the differential signal domain used to learn the probe-specific variance parameters. See model description in vignette for details.

tau2 A list. Each element corresponds to a probeset, and contains a vector that gives the estimated variance (noise level) for each probe in that probeset.

affinity A list. Each element corresponds to a probeset, and contains a vector that gives the estimated affinity effect for each probe in that probeset.

data Original probe-level data.

alpha, beta Inverse Gamma prior parameters for variance (tau2).

Extends

```
Class "list", from data part. Class "vector", by class "list", distance 2.
```

Author(s)

Leo Lahti <leo.lahti@iki.fi>

References

Probabilistic Analysis of Probe Reliability in Differential Gene Expression Studies with Short Oligonucleotide Arrays. Lahti et al., TCBB/IEEE 2011. http://www.roihu.info/publications/preprints/TCBB10.pdf

Examples

```
showClass("rpa.fit")
```

RPA.iteration

Estimating model parameters d and tau2.

Description

Finds point estimates of the model parameters d (estimated true signal underlying probe-level observations), and tau2 (probe-specific variances).

Usage

```
RPA.iteration(S, epsilon = 1e-3,
alpha = NULL, beta = NULL,
tau2.method = "fast", d.method = "fast",
maxloop = 1e6)
```

18 RPA.iteration

Arguments

S Matrix of probe-level observations for a single probeset: samples x probes.

epsilon Convergence tolerance. The iteration is deemed converged when the change in

all parameters is < epsilon.

alpha, beta Priors for inverse Gamma distribution of probe-specific variances. Noninforma-

tive prior is obtained with alpha, beta -> 0. Not used with tau2.method 'var'. Scalar alpha and beta are specify equal inverse Gamma prior for all probes to

regularize the solution. The defaults depend on the method.

tau2.method Optimization method for tau2 (probe-specific variances).

"robust": (default) update tau2 by posterior mean, regularized by informative priors that are identical for all probes (user-specified by setting scalar values for alpha, beta). This regularizes the solution, and avoids overfitting where a single probe obtains infinite reliability. This is a potential problem in the other tau2 update methods with non-informative variance priors. The default values alpha

= 2; beta = 1 are used if alpha and beta are not specified.

"mode": update tau2 with posterior mean "mean": update tau2 with posterior mean

"var": update tau2 with variance around d. Applies the fact that tau2 cost func-

tion converges to variance with large sample sizes.

d.method Method to optimize d.

"fast": (default) weighted mean over the probes, weighted by probe variances

The solution converges to this with large sample size.

"basic": optimization scheme to find a mode used in Lahti et al. TCBB/IEEE;

relatively slow; this is the preferred method with small sample sizes.

maxloop Maximum number of iterations in the estimation process.

Details

Assuming data set S with P observations of signal d with Gaussian noise that is specific for each observation (specified by a vector tau2 of length P), this method gives a point estimate of d and tau2. Probe-level variance priors alpha, beta can be used with tau2.methods 'robust', 'mode', and 'mean'. The d.method = "fast" is the recommended method for point computing point estimates with large sample size.

Value

A list with the following elements:

d A vector. Estimated 'true' signal underlying the noisy probe-level observations.

tau2 A vector. Estimated variances for each measurement (or probe).

Author(s)

Leo Lahti <leo.lahti@iki.fi>

References

Probabilistic Analysis of Probe Reliability in Differential Gene Expression Studies with Short Oligonucleotide Arrays. Lahti et al., TCBB/IEEE 2011.

rpa.list-class 19

Examples

```
## Not run:

## Preprocess probe-level data
## cind determines the 'reference' array
#Smat <- RPA.preprocess(Dilution, cind = 1)

## Pick probe-level data for one probe set
#pmindices <- pmindex(Dilution, "1000_at")[[1]]
#S <- t(Smat$fcmat[pmindices, ])

## RPA with default parameters:
#res <- RPA.iteration(S)
```

rpa.list-class

Class "rpa.list"

Description

Class for the RPA package.

Objects from the Class

Objects can be created by calls of the form new("rpa.list", ...).

Slots

An extended list. Contains the following information for one probeset.

Object of class "list"

- .DdtaA vector (probeset x arrays). Specifies the estimated 'true' underlying differential gene expression signal over the arrays (vs. the reference array 'cind') for the investigated probeset. Note that the reference array is not included.
- **mu.real** Deviation between the signal shape and the real signal level. Here d is the sum of mu.real and shape; affinities are calculated around d = mu.real + shape. Note that in Lahti11 paper d denotes the shape; we will polish the notation later.
- tau2 Contains a vector that gives the estimated variance for each probe in the investigated probeset.
- **affinity** A list. Each element corresponds to a probeset, and contains a vector that gives the estimated affinity effect for each probe in that probeset.
- **cind** Specifies which of the arrays in abatch was used as the reference for computing probe-level differential expression.
- set Probeset name.
- data Preprocessed probe-level data on which the model was fitted.

Extends

```
Class "list", from data part. Class "vector", by class "list", distance 2.
```

20 rpa.online

Methods

```
plot signature(x = "rpa.list"): ...
```

Author(s)

Leo Lahti < leo.lahti@iki.fi>

References

Probabilistic Analysis of Probe Reliability in Differential Gene Expression Studies with Short Oligonucleotide Arrays. Lahti et al., TCBB/IEEE 2011.

Examples

```
showClass("rpa.list")
```

rpa.online

RPA-online for preprocessing very large expression data sets.

Description

rpa.online is used to preprocess very large expression data collections based on a Bayesian hyperparameter update procedure. Returns an expressionSet object preprocessed with RPA.

Usage

```
\label{eq:colline} \begin{split} & \text{rpa.online}(\text{cel.path} = \text{NULL, cel.files} = \text{NULL, sets} = \text{NULL, cdf} = \\ & \text{NULL, bg.method} = \text{"rma", priors} = \text{list}(\text{alpha} = 2, \text{beta} = 1), \text{epsilon} = \text{1e-2, mc.cores} = 1, \\ & \text{verbose} = \text{TRUE, shuffle} = \text{TRUE, batch.size} = 100, \\ & \text{batches} = \text{NULL, quantile.basis} = \text{NULL, save.batches} = \text{TRUE,} \\ & \text{save.batches.dir, keep.batch.files} = \text{FALSE, unique.run.identifier} = \text{NULL}) \end{split}
```

Arguments

cel.path

cel.files	List of CEL files to preprocess
sets	Probesets for which RPA will be computed
cdf	Specify an alternative CDF environment
bg.method	Specify background correction method. See

Path to CEL file directory

bg.method Specify background correction method. See bgcorrect.methods() for options.

priors An 'rpa.priors' object. Can be used to set user-specified priors for the model

parameters. Not used tau2.method = "var". The prior parameters alpha and beta are prior parameters for inverse Gamma distribution of probe-specific variances. Noninformative prior is obtained with alpha, beta -> 0. Not used with tau2.method 'var'. Scalar alpha and beta specify an identical inverse Gamma prior for all probes, which regularizes the solution. Can be also specified as

lists, each element corresponding to one probeset.

epsilon Convergence tolerance. The iteration is deemed converged when the change in

all parameters is < epsilon.

mc.cores Number of cores for parallel computation

rpa.online 21

verbose Print progress information during computation

shuffle Form random batches

batch.size Batch size for online mode (rpa.online); the complete list of CEL files will be

preprocessed in batches with this size using Bayesian online-updates for probe-

specific parameters.

batches User-defined CEL file batches

quantile.basis Pre-calculated basis for quantile normalization

save.batches Save batches to speed up preprocessing.
save.batches.dir Output directory for temporary batch saves.

keep.batch.files Logical. Keep (TRUE) or remove (FALSE) the batch files after preprocessing.

unique.run.identifier

Define identifier for this run for naming the temporary batch files. By default, a

random id is generated.

Details

RPA preprocessing function. Gives an estimate of the probeset-level mean parameter d of the RPA model, and returns these in an expressionSet object. The CEL files are handled in batches to obtain Bayesian updates for probe-specific hyperpriors; after sweeping through the database in batches the results are combined. The online mode is useful for preprocessing very large expression data sets where ordinary preprocessing algorithms fail, without compromises in modelling stage.

Value

An instance of the 'expressionSet' class.

Note

Development version

Author(s)

Leo Lahti leo.lahti@iki.fi>

References

See citation("RPA").

See Also

rpa, AffyBatch, ExpressionSet

Examples

```
# Not run:
# eset <- rpa.online(cel.file.path)
```

22 rpa.plot

rpa.plot	Plot RPA results and probe-level data for a specified probeset.

Description

Plots the preprocessed probe-level observations, estimated probeset-level signal, and probe-specific variances. It is also possible to highlight individual probes and external summary measures.

Usage

Arguments

dat Original data: probes x samples. rpa.fit.object An instance of the 'rpa.fit' class.

toydata.object Optional. Output from sample.probeset toydata generator function. Can be

used to compare (toy)data with known ground truth to RPA estimates from

rpa.fit.object.

highlight.probes

Optionally highlight some of the probes (with dashed line)

pcol Color for probe signal visualization.

mucol Color for summary estimate.
ecol Color for external signal.
cex.lab, cex.axis, cex.main, cex.names

Font size adjustment parameters.

external.signal Plot external signal on the probeset. For instance, an alternative summary esti-

mate from another preprocessing methods

main Title text.

plots "all": plot data and summary, noise and affinity "data": plot data and summary

... Other parameters to pass for plot function.

Value

Used for its side-effects. Returns probes x samples matrix of probe-level data plotted on the image.

Author(s)

Leo Lahti@iki.fi>

References

Probabilistic Analysis of Probe Reliability in Differential Gene Expression Studies with Short Oligonucleotide Arrays. Lahti et al., TCBB/IEEE. See http://www.cis.hut.fi/projects/mi/software/RPA/

RPA.pointestimate 23

See Also

RPA.pointestimate

Examples

```
# Not run:

## Load example data set
#require(affydata)
#data(Dilution)

## Compute RPA for specific probesets only
#set <- "1000_at"
#rpa.results <- RPA.pointestimate(Dilution, set)

## Visualize the results for one of the probe sets
#rpa.plot(set, rpa.results)
```

RPA.pointestimate

Computing point estimate for the model parameters for all probe sets.

Description

Computes point estimate

Usage

```
RPA.pointestimate(abatch, sets = NULL, myseed = 101, priors = NULL, epsilon = 1e-2, cind = 1, tau2.method = "robust", d.method = "fast", verbose = TRUE, bg.method = "rma", normalization.method = "quantiles.robust", cdf = NULL)
```

Arguments

abatch An AffyBatch object.

sets Specifies the probesets for which RPA estimates will be computed. Default: all

probe sets.

myseed Specifies the random seed.

priors Optional list containing hyperparameters alpha and beta of the inverse Gamma

prior of the probe-specific variances; alpha is a scalar, common for all probes and probeset; beta is a list where each element is a vector corresponding to one probeset, specifying beta for each probe. Can be used to set user-specified priors for the model parameters. Not applicable for tau2.method = "var". Noninformative prior is obtained with alpha, beta -> 0. Not used with tau2.method 'var'. Can be used to regularize the solution with small sample size, or in batch-wise online-updates with large sample size. If priors are not provided for certain

probesets (NULL), default priors are used.

epsilon Convergence tolerance. The iteration is deemed converged when the change in

the d parameter is < epsilon.

24 RPA.pointestimate

cind Specifies which array in abatch is used as a reference in computing probe-level

differential expression.

tau2.method Optimization method for tau2 (probe-specific variances).

"robust": (default) update tau2 by posterior mean, regularized by informative priors that are identical for all probes (user-specified by setting scalar values for alpha, beta). This regularizes the solution and avoids overfitting where a single probe obtains infinite reliability. This is a potential problem in the other tau2 update methods with non-informative variance priors. The default values alpha

= 2; beta = 1 are used if alpha and beta are not specified.

"mode": update tau2 with posterior mean "mean": update tau2 with posterior mean

"var": update tau2 with variance around d. Applies the fact that tau2 cost func-

tion converges to variance with large sample sizes.

d.method Method to optimize d.

"fast": (default) weighted mean over the probes, weighted by probe variances

The solution converges to this with large sample size.

"basic": optimization scheme to find a mode used in Lahti et al. TCBB/IEEE;

relatively slow; this is the preferred method with small sample sizes.

verbose Print progress information during computation. Default: TRUE.

bg.method Specify background correction method. Default: "rma". See bgcorrect.methods()

for other options.

normalization.method

Specify quantile normalization method. Default: "pmonly". See normalize.methods(Dilution)

for other options.

cdf Specify an alternative CDF environment. Default: none.

Details

Calculates RPA estimates of probe reliability and differential expression between the user-specified reference array (cind) and the other arrays in the data set. The model assumes P observations for each transcript target (i.e. a probeset) with Gaussian noise which is specific for each probe (variance is specified by tau2). The mean (affinity) parameters of the Gaussian noise model cancel out in calculating probe-level differential expression. RPA.pointestimate gives a point estimate for d and tau2. The 'prior' parameter is not applicable with tau2.method = "var". The d.method = "fast" is recommended with large sample size.

Value

An instance of class 'rpa'. This is an extended list containing the following elements:

d A matrix of probesets x arrays. Specifies the estimated 'true' underlying differ-

ential gene expression signal over the arrays (vs. the reference array 'cind') for

each investigated probeset. Note that the reference array is not included.

tau2 A list. Each element corresponds to a probeset, and contains a vector that gives

the estimated variance for each probe in that probeset. This corresponds to the

parameter tau^2 in the vignette and manuscript.

cind Specifies which of the arrays in the abatch (the affybatch object to be analyzed)

has been used as the reference for computing probe-level differential expression.

affinity Probe affinity effects.

sets A character vector listing the investigated probesets.

RPA.preprocess 25

Note

tau2.method = "robust" and d.method = "fast" are recommended. With small sample size and informative priors, d.method = "basic" may be preferable, with large sample size d.method = "fast" should have considerable speedups and comparable accuracy.

Author(s)

```
Leo Lahti <leo.lahti@iki.fi>
```

References

```
See citation("RPA").
```

See Also

rpa.plot, rpa, set.priors, rpa2eset, RPA.preprocess, AffyBatch, rpa.fit, estimate.affinities

Examples

```
## Load example data set
#require(affydata)
#data(Dilution)

## Compute RPA for whole data set
## ... slow, not executed here
# rpa.results <- RPA.pointestimate(Dilution)

## Visualize the results for one of the probe sets
#rpa.plot(set, rpa.results)
```

RPA.preprocess

Preprocess AffyBatch object for RPA.

Description

Background correction, quantile normalization and log2-transformation for probe-level raw data in abatch. Then probe-level differential expression is computed between the specified 'reference' array (cind) and the other arrays.

Usage

```
\label{eq:RPA.preprocess} RPA.preprocess (abatch, bg.method = "rma", \\ normalization.method = "quantiles.robust", \\ cdf = NULL, \\ cel.files = NULL, cel.path = NULL)
```

26 rpa2eset

Arguments

abatch An AffyBatch object.

bg.method Specify background correction method. See bgcorrect.methods(abatch) for op-

tions.

normalization.method

Specify normalization method. See normalize.methods(abatch) for options. For

memory-efficient online version, use "quantiles.online".

cdf The CDF environment used in the analysis.

cel.files List of CEL files to preprocess.
cel.path Path to CEL file directory.

Details

Probe-specific variance estimates are robust against the choice of reference array.

Value

fcmat Probes x arrays preprocessed differential expression matrix.

cind Specifies which array in abatch was selected as a reference in calculating probe-

level differential expression.

cdf The CDF environment used in the analysis.

set.inds Indices for probes in each probeset, corresponding to the rows of fcmat.

Author(s)

Leo Lahti <leo.lahti@iki.fi>

References

See citation("RPA").

See Also

AffyBatch

rpa2eset Coerce 'rpa' object into an 'ExpressionSet'

Description

An instance of 'rpa' class contains differential gene expression estimates in the variable 'd'. The function 'rpa2eset' coerces this into an ExpressionSet object to allow downstream analysis of the results using standard R/BioC tools for gene expression data.

Usage

rpa2eset(x)

s2.update 27

Arguments

X

An instance of the rpa class (obtained as output from RPA.pointestimate)

Value

An 'ExpressionSet' object.

Author(s)

```
Leo Lahti <leo.lahti@iki.fi>
```

References

Probabilistic Analysis of Probe Reliability in Differential Gene Expression Studies with Short Oligonucleotide Arrays. Lahti et al., TCBB/IEEE. See http://bioconductor.org/packages/release/bioc/html/RPA.html

Examples

```
# Not run:

#require(RPA)
#require(affy)
#require(affydata)
#data(Dilution)

## Compute RPA for specific probesets
#sets <- geneNames(Dilution)[1:2]
#rpa.results <- RPA.pointestimate(Dilution,sets)

## Coerce the rpa object into an ExpressionSet
#eset <- rpa2eset(rpa.results)
```

s2.update

Variance update.

Description

Update variance parameters in the Bayesian RPA model with an EM-type optimization procedure and potentially informative priors.

Usage

```
s2.update(dat, alpha = 1e-2, beta = 1e-2, s2.init = NULL, th = 0.01)
```

Arguments

dat A probes x samples matrix (probeset).

alpha, beta Shape and scale parameters of the inverse gamma prior for the variances.

s2.init Initial values for the variances for optimization.

th Optimization convergence threshold.

28 sample.probeset

Details

Updates the variance parameters in the Bayesian RPA model with a (generalized) EM-type optimization procedure and potentially informative priors. The variances are updated by their mode at each step until convergence. Priors (alpha, beta) are taken into account if provided.

Value

A vector of variance parameters, one for each probe.

Author(s)

```
Leo Lahti <leo.lahti@iki.fi>
```

References

```
See citation("RPA").
```

See Also

rpa.online

Examples

```
\# dat: probes x samples matrix (probeset) \# s2 <- s2.update(dat)
```

sample.probeset

Toydata generator for probeset data.

Description

Generate random probeset with varying probe-specific affinities and variances.

Usage

```
sample.probeset(P = 10, n = 20, shape = 1, scale = 1, mu.real = 2)
```

Arguments

P Number of probes.n Number of samples.

shape, scale Shape and scale parameters of the inverse Gamma function used to generate the

probe-specific variances.

mu.real Absolute signal level of the probeset.

Details

The toy data generator follows distributional assumptions of the RPA model and allows quantitative estimation of model accuracy with different options, noise levels and sample sizes. Probeset-level summary estimate is obtained as mu-real + d.

summarize.batches 29

Value

A list with the following elements:

dat Probeset data: probes x samples

variance Probe variances. affinity Probe affinities.

d Probeset-level signal shape.

mu.real Probeset-level absolute signal level.

Author(s)

Leo Lahti <leo.lahti@iki.fi>

References

See citation("RPA").

Examples

```
real <- sample.probeset(P = 10, n = 20, shape = 1, scale = 1, mu.real = 2)
```

summarize.batches Summarize batches

Description

Probeset summarization in batches

Usage

```
summarize.batches (sets = NULL, variances, batches, load.batches = FALSE,\\ mc.cores = 1, cdf = NULL, bg.method = "rma", normalization.method = "quantiles", verbose = TRUE, quantile.basis, save.batches.dir = ".", unique.run.identifier = NULL, save.batches.
```

Arguments

sets Probesets to summarize

variances Precalculated probe-specific variances batches Data batches for online learning

load.batches Logical. Load precalculated data for the batches.

mc.cores Number of cores for parallel computation cdf CDF for alternative probeset definitions

bg.method Background correction method

 ${\bf normalization.} \\ {\bf method}$

Normalization method

verbose Print progress information

quantile.basis Basis for quantile normalization

save.batches.dir Specify the output directory for temporary batch saves.

30 summarize.batches

unique.run.identifier

Define identifier for this run for naming the temporary batch files. By default, a random id is generated.

save.batches Save batches?

Details

Sweeps through the batches. Summarizes the probesets within each batch based on the precalculated model parameters using point estimates of the model.

Value

Expression matrix: probesets x samples.

Author(s)

Leo Lahti <leo.lahti@iki.fi>

References

See citation("RPA").

Examples

#

Index

*Topic classes rpa-class, 14	plot,rpa.list-method (rpa.list-class), 19 probe.performance, 11
rpa-class, 14	probe.performance, 11
rpa.list-class, 19	qnorm.basis.online (online.quantile), 9
*Topic iteration	quorimisasisionimis (ommisiquanono), y
-	RPA (RPA-package), 2
hyperparameter.update, 7	rpa, 12
RPA.iteration, 17	rpa-class, 14
s2.update, 27	RPA-package, 2
*Topic methods	rpa.fit, 15, 16
get.probe.noise.estimates, 6	-
plot, 10	rpa.fit-class, 16
rpa, 12	RPA.iteration, 17
RPA.iteration, 17	rpa.list-class, 19
rpa.online, 20	rpa.online, 20
rpa.plot, 22	rpa.plot, 22
RPA.pointestimate, 23	RPA.pointestimate, 14, 23
RPA.preprocess, 25	RPA.preprocess, 25
*Topic package	rpa2eset, 26
RPA-package, 2	
*Topic utilities	s2.update, 27
estimate.affinities, 3	sample.probeset, 28
estimate.hyperparameters, 4	show,rpa-method (rpa-class), 14
get.batches, 5	summarize.batches, 29
online.quantile, 9	
	vector, 14, 17, 19
probe.performance, 11	
rpa.fit, 15	
rpa2eset, 26	
sample.probeset, 28	
summarize.batches, 29	
[,rpa-method (rpa-class), 14	
[[,rpa-method (rpa-class), 14	
estimate.affinities, 3	
estimate.hyperparameters, 4	
get.batches, 5	
get.probe.noise.estimates, 6	
hyperparameter.update, 7	
list, 14, 17, 19	
online.quantile, 9	
ommo. quantito, y	
plot, 10	